

Refraction unit MINI Master



MINI Master



Refraction unit MINI Master

- Sliding table top for two instruments.
- Integrated adjustable power transformer for slit lamp/ophthalmometer.
- Adapter for chin rest.
- Table tops made using laminated wood.
- To be fixed to the floor.

Optional

- Left-handed version.
- Factory installed data cables.
- Foot switch control for up/down movement of the table.



Facilities for
the disabled

Refraction unit

Med
instrus

Refraction unit MINI Master



1. Adjustable power supply for instruments.

2. Sliding table top for two instruments.



3. Control panel.

Technical information

Power supply	220-230V/50Hz
Max weight	80 kg
Table height from the ground	From 70 cm to 95 cm (lift 25 cm)
Adjustable power supply for instruments	6V and/or 12V (other voltage on request)

