


# PORTABLE SLIT LAMP LS-1A/B





## ◆ EXCELLENT OPTICAL QUALITY

Especially fit for children, ward round, screening

## ◆ BIG SLIT ANGLE

Illumination angle max. up to be 60°, can get better slit view for cornea etc.

## ◆ WIDE VISUAL FIELD

50% higher field of view than traditional design and a wider observation range

## ◆ ONE-HAND OPERATION

Easily adjust observation magnification during the inspection, help doctors quickly find the lesion in a larger field of vision and transfer to a higher magnification for better observation

## ◆ WARM LED

Warm light for better clinical observation of eye tissues and LED light source for longer and safer use

## ◆ SPECIFICATION

|                      |                                            |
|----------------------|--------------------------------------------|
| Magnification        | LS-1A:10X<br>LS-1B:10X/16X                 |
| Working Distance     | 80mm(16X)<br>100mm(10X)                    |
| Field of View        | ø10mm(16X)<br>ø15mm(10X)                   |
| Interpupillary       | 49-75mm                                    |
| Diopter Adjustment   | ±7D                                        |
| Slit Size            | 0.1mm*12mm<br>0.2mm*8mm<br>0.8mm*12mm      |
| Spot Size            | ø1 / ø5 / ø12mm                            |
| Slit Angel           | Horizontal, ±60°                           |
| Filter               | Neutral Density<br>Red-free<br>Cobalt Blue |
| Light Resource       | Warm LED                                   |
| Illumination control | Continuously<br>Adjustable                 |
| Battery Type         | 2 AA batteries                             |
| Operating Time       | ≥4 Hours                                   |
| Dimensions           | 320*310*205mm                              |
| Net Weight           | Approx. 835g                               |



**ChongQing Sunkingdom Medical Instrument Co., Ltd.**

Head office: Room 12, Floor 10, Block A of China Resource Center. Jiulongpo District, Chongqing 400050, China.  
 FAX: +86 23 68102793  
 Tel: +86 23 68643990

E-mail: [sales@cqsunkingdom.com](mailto:sales@cqsunkingdom.com)  
 URL: [en.cqsunkingdom.com](http://en.cqsunkingdom.com)