


Cyclops Lite

The easiest to use iPad dispensing measuring device ever!



RPD/LPD	PD	BovA	BovB	REP/LEP	BovED	DBL
31.0/31.2	62.2	52.2	31.5	16.3/14.8	60.9	16.5
VCD	PT	ANPD/LNPD	NPD	PD		
8.0	3.1	25.3/20.6	54.9	33		

Plastic Frame



- ▶ Distance PD measurement
- ▶ Fitting height/Wrap angle/PT/VCD measurement
- ▶ Near PD & Reading distance measurement
- ▶ Boxing data

CYCLOPS Lite

Key feature

- ▶ Distance PD and Near PD & Reading distance measurement
- ▶ Fitting height/Wrap angle/PT/VCD measurement
- ▶ Flash light/Auto light adjustment
- ▶ Wide range of measuring distance 50cm~64cm
- ▶ Quick & Easy to attach and detach
- ▶ Extremely light (Weight only 148g)
- ▶ Optional Desk stand for ease of use
- ▶ Convex lens for patient fixation
- ▶ Accurate/Repeatable results


Measuring feature

Distance PD measurement	Fitting height/Wrap angle/PT/VCD
Near PD & Reading distance measurement	Boxing data

Other feature

Flash light/Auto light adjustment	Manual & Auto capture
Accurate & Repeatable results	Wide range of measurement distance 50cm~64cm
Ease of use by using optional Desk stand	Quick & Easy to install/remove
Extremely light (Weight only 148g)	Convex lens for patient fixation

Measurement procedure



DESK STAND

Auxiliary Equipment

The optional Desk stand helps maintain a stable posture during measurement. It is fully adjustable and quick to install & remove.

- Accurate results/Adjustable/Stability for shooting



Product specification

- Color : White • Mirror : Convex lens

Size and weight

Size : 180 x 88 x 20mm
Weight : 148g

Power and battery

Charging via USB to computer system or other outlet
Fully charged within 4 hours (if 300mA supplied)

System requirement

iPad IOS 6.0 or later
Cyclops_Lite application(will be downloadable from App store)
No jailbreak iPad

In the package

Cyclops_Lite /
USB charging cable /
Jig / "L" wrench

Eye-care business venture company

VIEW W I T E C H
www.viewitech.com

We have special price for optical shop franchiser/product alliance, Please contact us.
Tel : +82-31-8086-7650 Fax : +82-31-8086-7660
E-mail: sales@viewitech.com
Copyright© Viewitech Co.,Ltd. All Rights Reserved